

3D Printing Technology for Clinical Applications



James J. Yoo

**Wake Forest Institute for Regenerative Medicine
Wake Forest School of Medicine
Winston-Salem, North Carolina**




Strategy for Clinical Translation


Cell



Scaffold




Patient



- ❖ **Cells alone – Utilize or Enhance Cellular Function**
 - Somatic Cells, Stem Cells
- ❖ **Scaffolds alone – Bridging tissue defects, Regeneration**
 - Natural, Synthetic, Hybrid Biomaterials
- ❖ **Cells + Scaffolds – Engineered Tissues for Regeneration**


Tissue Engineering of the Bladder



TRANSLATION

COMPLETED

Urethra



Vagina




BI Neck



Ureter



Bladder



CURRENT DEVELOPMENT

Blood Vessel




Heart Valve



Liver



Corpora




Ear




Finger




Kidney




Nerve



Skin




Muscle




- **Tissue engineering has had initial successes with building a number of tissues clinically**
- **Challenges still exist in developing and translating complex tissue systems that require coordinated and systemic function**

Engineering Challenge



Vascular Network




Organ and tissue complexity




Bioprinting Technology

- Bioprinting technology is evolving into a viable method to fabricate tissue constructs for reconstruction
- Bioprinting technology can deliver 1) multiple cell types, 2) gel biomaterials, and 3) genes, proteins and macromolecules to target sites
- Multi-cellular tissue/organ constructs can be generated simultaneously with the maintenance of their structural and spatial orientation *in vivo*

Layer-by-layer assembly
of 3-D constructs






Bioprinting Technology



Heart

Bone

Blood Vessels Composite Tissue






Skin Printer for Burn



- Burns comprise 5-20% of injuries in combat zones
- Large area burns are difficult to stabilize in the combat theater

Proposed Solution:

- Develop a novel delivery system that would repair battlefield burn injuries *in situ* using a portable skin printing system



Design and Construction of Clinical Skin Bioprinter



Supported, in part, by TATRC

Validation Study Design


Skin Biopsy





Cell Isolation and Expansion



100 cm² Full-Thickness Wound (n=6)



Cell Delivery



Suspension of 10⁶ cells in solution

MAT Only Cell Printing



Untreated Cell Spraying

Data Collection & Analyses




- Gross observation
 - Wound closure
 - Contracture
 - Re-epithelialization
- Histology
 - H & E
 - Masson Trichrome
 - Immunohistochemistry

Printing Process


Skin Wound




Scanning Wound




Raw Image




Standard Coordinates



Printed Wound




Bioprinting



Operating Software



Generation of Nozzle Path



Morphological Analysis





Untreated

Matrix Only





Cell Spray

Cell Printed

2 wks






4 wks




Wounds bioprinted with cells showed well organized re-epithelialization by 2 weeks

- **Generation of 3D freeform shaped constructs with precision**
 - Multiple cell types, biomaterials, drugs
- **High strength constructs:**
 - Gel and polymeric materials (~12)
- **Printing resolution:**
 - Cell printing: $\geq 50 \mu\text{m}$
 - Structural material printing: $\geq 2 \mu\text{m}$





Working Principle




Fabrication Process





- *White: Polycaprolactone*
- *Pink: Cell mixed Hydrogel*

3-D Cell Patterning


Red/Green
 Color Cells PCL Sacrificial
 material




Printing w/
 Synthetic polymer
 & fluorescent labeled
 cells



PCL
 Channel
 Red
 color cell
 Green
 color cell





Cell Viability




Fluorescence Images


NIH3T3 Cells
Live Cell (Green) / Dead Cell (Red)




Liver Structure




Bladder Tissue Structure




Testis









Cardiac Muscle






Kidney Structure




Fabrication of Engineered Constructs: Bone




Fabrication of Muscle Constructs



In vitro




In vivo




Institute for Regenerative Medicine *Chondrocyte laden Cell Carrier Material*




Ear CAD model




Scanning Path



Printing Condition








In Vitro



In Vivo




Muscle-Tendon Junction (MTJ)



Red: Desmin; Green: Collagen type I

Integrated Organ Printing System


Reverse Engineering




Medical Data




3D CAD Model



CAD/CAM System



Stage Operation Code



Fabrication System



3D Structure

Code Generation

Courtesy of Prof. DW Cho, Postech

Conclusion

- **Bioprinting technology is able to generate 3-D tissue constructs with precision**
- **Delivered cells are uniformly distributed within the construct and remain viable over time**
- **This system can be used as a cell delivery tool for building complex and composite tissue and organ systems**



- **Technology Development**
Hardware, Software, Bioink (cells, materials)
- **Regulatory Pathway**
- **Commercialization Strategy**

Wake Forest Institute for Regenerative Medicine



- **WFIRM: >300 Members: (Dr. Anthony Atala-Director)**
- **Working on over 30 different tissues and organs**

This work was made possible, in part, by grants from the following institutions:

NIH: NIDDK

NIH: HLI

Department of Defense (AFIRM, TATRC, OETRP, NHRC)

Department of Energy

National Kidney Foundation

Muscular Dystrophy Association

The Crown Foundation

The Frase Foundation

The Nakos Foundation

The Chapman Foundation

JDRF

Musculoskeletal Transplant Foundation

Tengion, Inc

Plureon

Errett Fisher Foundation